
**Tube set for insufflation, single-use / Schlauchset für Insufflation, einmalverwendbar / Комплект маркучи за инсуфлиране, еднократна употреба / Sada hadiček na insuflaci, na jedno použití / Slangesæt til insufflation, engangsbrug / Lot de tubulures pour insufflation, usage unique / Σετ εύκαμπτων σωλήνων για εμφύσηση, μίας χρήσης / Tömlőkészlet inszufflációhoz, egyszer használatos / Set di tubi per insufflazione, monouso / Insuflācijas caurules komplekts, vienreizējas lietošanas / Insuflācijas lanksčiujų vamzdelių komplektas, vienkartinis / Slangesett til insufflasjon, til engangsbruk / Zestaw drenów do insuflacji, jednorazowy / Jogo de tubos para insuflação, descartável / Set furtun pentru insuflație, de unică folosință / Hadicová súprava na insufláciu, jednorazová / Komplet cevi za insuflacijo za enkratno uporabo / Juego de tubos para insuflación, desechable / Slangsats för insufflering, för engangsbruk**

**CE 0197**

**W.O.M. WORLD OF MEDICINE GmbH**  
Salzufer 8  
10587 Berlin, Germany  
Phone: +49 30 39981-550  
Fax: +49 30 39981-545  
E-mail: info.berlin@novanta.com

CE marking according to Regulation (EU) 2017/745  
CE-Kennzeichnung gemäß Verordnung (EU) 2017/745  
CE маркировка съгласно Регламент (ЕС) 2017/745  
Značka CE podle nařízení (EU) 2017/745  
CE-mærkning iht. forordningen (EU) 2017/745  
Marquage CE selon le règlement (UE) 2017/745  
Σήμανση CE σύμφωνα με τον Κανονισμό (ΕΕ) 2017/745  
CE jelölés a 2017/745EU rendelet szerint  
Marcatura CE ai sensi del Regolamento (UE) 2017/745  
CE marķējums atbilstoši Regulai (ES) 2017/745  
CE ženklas pagal Reglamentą (ES) 2017/745  
CE-merking i samsvar med direktiv (EU) 2017/745  
Oznaczenie CE zgodnie z Rozporządzeniem (UE) 2017/745  
Marcação CE de acordo com o regulamento (UE) 2017/745  
Marcaj CE în conformitate cu Regulamentul (UE) 2017/745  
CE označenie podľa smernice (EÚ) 2017/745  
Oznaka CE v skladu z Direktivo (EU) 2017/745  
Marcado CE conforme al Reglamento (UE) 2017/745  
CE-märkning enligt förordning (EU) 2017/745

**Instructions for Use**

# ST293

en

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# ST293

de

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sv

## Symbols and Descriptions

|                                                                                     |                                                                |                                                                                      |                                           |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------|--------------------------------------------------------------------------------------|-------------------------------------------|
|    | Consult instructions for use                                   |   | Use by date (YYYY-MM-DD)                  |
|    | Medical Device                                                 |    | Quantity                                  |
|    | Model number                                                   |    | Do not use if package is damaged          |
|    | Unique device identifier                                       |    | Keep dry                                  |
|    | Authorized for Sale or use by Physician only                   |    | Keep away from sunlight                   |
|    | Caution                                                        |   | Protect from heat and radioactive sources |
|   | Catalogue number                                               |  | Humidity limitation                       |
|  | Batch code                                                     |  | Atmospheric pressure limitation           |
|  | Manufacturer                                                   |  | Temperature limit                         |
|  | Country (DE) and date of manufacture (YYYY-MM-DD)              |  | Storage conditions                        |
|  | Single sterile barrier system, sterilized using ethylene oxide |  | Transport conditions                      |
|  | Sterilized using ethylene oxide                                |  | Not made with phthalates                  |
|  | Do not resterilize                                             |                                                                                      |                                           |
|  | Do not reuse                                                   |                                                                                      |                                           |

## 1 Important User Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing or service personnel, or
- damages or malfunction of device and/or accessories.

To request further copies of the instructions for use, please contact the manufacturer.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

Paragraphs marked with the words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be given special attention.

**Subject to technical changes**

**Safety and warning messages**

### WARNING!

The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third party.

### CAUTION!

These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

### NOTE!

These paragraphs contain information to clarify the instructions or provide additional useful information.


## 2 Description

The tube set is an accessory part of the FM134 CO<sub>2</sub> insufflator for endoscopic procedures. The sterile tube set serves to transport insufflation gas from the FM134 CO<sub>2</sub> insufflator to the patient.

The tube set is used with the FM134 CO<sub>2</sub> insufflator with an ISO connector.

**Fig. 2-1 Parts of the tube set**

- ① ISO connector
- ② Filter
- ③ Tube line
- ④ Luer lock connector
- ⑤ Protection cap

**Target patient group** The tube set is used with patients for whom a laparoscopy or an endoscopy vessel harvesting procedure is indicated. The target patient group is not limited by age, gender, or condition of health. The tube set may not be used with patients weighing less than 2.5 kg.

**2.1 Intended Use**  
The tube set for insufflation, single-use is intended to transport CO<sub>2</sub> gas from the insufflator to the patient.

**2.2 Indications and Contraindications**  
**Indications** Indications are identical to the FM134 CO<sub>2</sub> insufflator. Please refer to the Instructions for Use of the FM134 CO<sub>2</sub> insufflator.

**Contraindications** Contraindications are identical to the FM134 CO<sub>2</sub> insufflator. Please refer to the Instructions for Use of the FM134 CO<sub>2</sub> insufflator.

**2.3 General Warnings CO<sub>2</sub> Insufflation**  
Regarding the residual risks related to CO<sub>2</sub> insufflation, please refer to the Instructions for Use of the FM134 CO<sub>2</sub> insufflator:

- Idiosyncratic reactions
- CO<sub>2</sub> absorption
- Metabolic and cardiac reactions/ Metabolic, circulatory and respiratory reactions
- Hypothermia/monitoring body temperature
- Dehydration
- Embolism/insufflation of internal organs
- Subcutaneous emphysema

**3 Setup and Usage**

**WARNING!**  
**Check for defects and expiration date**  
Check the product and the packaging for defects prior to use. Make sure that the expiration date has not been exceeded.  
Do not use if package is damaged or if the expiration date has been exceeded.

**WARNING!**  
**Leaky insufflation tube**  
Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.

**WARNING!**  
**Contamination**  
Do not use the tube set if signs or indications of a contamination are detected.

**WARNING!**  
**Use by healthcare professional staff only**  
The device must only be used by surgeons or operating room personnel with the necessary professional qualifications.

**WARNING!****Professional qualification**


The instructions for use do not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

**CAUTION!****In case of serious incident**

Please report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the State in which the user and/or patient is established.

Fig. 3–1 Connecting tube set

- ① FM134 CO<sub>2</sub> insufflator
- ② ISO connector, FM134 CO<sub>2</sub> insufflator
- ③ ISO connector, tube set
- ④ Tube set with integrated filter

**ISO connection**

An ISO connector is located at the FM134 CO<sub>2</sub> insufflator. Connect the ISO connector of the tube set firmly with the ISO connector of the FM134 CO<sub>2</sub> insufflator.

## 4 Connecting the Tube Set

**WARNING!****Air in the insufflation tube set**

Air in the insufflation tube set can cause an embolism.

Start the insufflation before the tube set is connected to the patient to push the air out of the insufflation tube set. Stop the insufflation when at least 1 liter of CO<sub>2</sub> has been insufflated.

**NOTE!****Limitation of inflow**

Kinking or twisting of the tubes may lead to limitation of inflow or to loosening of the Luer lock connection.

When placing tube sets on the OR table or the patient, make sure to allow for sufficient movement of the tubes.

When connecting the Luer lock to a trocar, make sure it is secured tightly.

Removing the tube set from the packaging

Connect Veress needle or trocar

Connect tube set to CO<sub>2</sub> insufflator

Function test of the tube set

**1. To be carried out by non-sterile personnel:**

- Open packaging of the tube set and deliver the tube set to the sterile personnel.

**2. To be carried out by sterile personnel:**

- Keep the Luer lock connector in the sterile area and hand the tube end with the filter to the non-sterile personnel.
- Remove the protection cap from the Luer lock connector.
- Connect the Luer lock connector with the Veress needle or trocar. Close the stopcock at the Veress needle or trocar.

**3. To be carried out by non-sterile personnel:**

- Turn the FM134 CO<sub>2</sub> insufflator on.
- Connect the ISO connector of the tube set to the ISO connector of the FM134 CO<sub>2</sub> insufflator.

**4. To be carried out by non-sterile and sterile personnel:**

**Non-sterile personnel:**

- Set the nominal pressure to 15 mmHg and the gas flow to 3 l/min. Start insufflation.
- Check the gas consumption on the display of the FM134 CO<sub>2</sub> insufflator. If the gas consumption is more than approximately 0.4 l after 30 seconds a leakage is indicated.

**Sterile and non-sterile personnel:**

- If a leakage is detected, repeat without Veress needle or trocar, manually sealing the tube end. If no leakage is now present, Veress needle or trocar must be replaced.
- If a leakage is still detected, repeat without tube set, manually sealing the FM134 CO<sub>2</sub> insufflator. If no leakage is now present, the tube set must be replaced.
- If a leakage is still detected, the FM134 CO<sub>2</sub> insufflator may be faulty, please refer to the Instructions for Use of the FM134 CO<sub>2</sub> insufflator.
- Ensure no leak can be detected before starting the operation.
- Allow at least 1 l of CO<sub>2</sub> to flow through the tube set to remove the air from the tube set before connecting the Veress needle or trocar to the patient.
- The tube set is now ready for use. Disconnect the tube set from the Veress needle or trocar before the Veress needle or trocar is inserted into the patient. The tube set can then be reconnected to the Veress needle or trocar and the insufflation started.

**5 Disconnection and Disposal of the Tube Set**

**WARNING!**

**Reprocessing of sterile disposable products**

Reuse of sterile disposable products can lead to patient and/or user infection hazards. Product functionality might be impaired. Risk of injury, illness, or death due to contamination and/or impaired functionality of the product. Do not reprocess the product.

**NOTE!**

**Disposing tube sets**

Observe applicable hygiene rules and regulations when disposing tube sets, collected fluids, and waste containers.

Remove the tube set only after insufflation is stopped. See the steps for stopping insufflation in the instructions for use of the FM134 CO<sub>2</sub> insufflator.

## 6 Technical Data

| General information                         |                                   |                                                                                                                                                          |
|---------------------------------------------|-----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product name                                |                                   | Tube set for insufflation, single-use                                                                                                                    |
| Model                                       |                                   | ST293                                                                                                                                                    |
| REF/catalogue number                        |                                   | Z0293                                                                                                                                                    |
| General product specifications              |                                   |                                                                                                                                                          |
|                                             | Sterilization                     | EO sterilization                                                                                                                                         |
|                                             | Weight                            | 150 g                                                                                                                                                    |
|                                             | Shelf life                        | 3 years                                                                                                                                                  |
|                                             | Flow                              | Min. 15 l/min                                                                                                                                            |
| Tubing                                      |                                   |                                                                                                                                                          |
| Tube line                                   | Total length                      | 3 m ± 0.1                                                                                                                                                |
|                                             | Tube diameter                     | 8.8 mm ± 0.15                                                                                                                                            |
|                                             | Material                          | PVC                                                                                                                                                      |
|                                             | Color                             | Transparent                                                                                                                                              |
| Connector (proximal)                        | Type                              | ISO connector 15 mm, male                                                                                                                                |
|                                             | Material                          | MABS                                                                                                                                                     |
|                                             | Color                             | Grey                                                                                                                                                     |
| Connector (distal)                          | Type                              | Luer lock connector, male                                                                                                                                |
|                                             | Material                          | MABS                                                                                                                                                     |
|                                             | Color                             | Transparent/blue                                                                                                                                         |
| Protection cap                              | Type                              | Vented cap for Luer lock                                                                                                                                 |
|                                             | Material                          | PE                                                                                                                                                       |
|                                             | Color                             | Transparent                                                                                                                                              |
| Functional components                       |                                   |                                                                                                                                                          |
| Filter                                      | Filter effectivity of filtermedia | BFE ≥ 99.9999 %<br>VFE ≥ 99.9999 %                                                                                                                       |
|                                             | Material                          | MABS (housing)<br>Hydrophobic glass-fibre media (filter media)                                                                                           |
|                                             | Color                             | Transparent/grey (housing)<br>White (filter media)                                                                                                       |
| Packaging                                   |                                   |                                                                                                                                                          |
| Sterile packaging                           |                                   | Single Tyvek blister (295 x 315 mm)                                                                                                                      |
| Box                                         |                                   | 20 blisters (580 x 240 x 310 mm)                                                                                                                         |
| Transport, storage and operating conditions |                                   |                                                                                                                                                          |
| Transport conditions                        |                                   | -18 to 55 °C / 0 to 131 °F<br>10 to 85 % rel. humidity, non-condensing<br>70 to 106 kPa air pressure                                                     |
| Storage conditions                          |                                   | 0 to 35 °C / 32 to 95 °F<br>10 to 85 % rel. humidity, non-condensing<br>70 to 106 kPa air pressure                                                       |
| Operating conditions                        |                                   | 10 to 30 °C / 50 to 86 °F<br>30 to 75 % rel. humidity, non-condensing<br>70 to 106 kPa air pressure<br>3000 m / 9843 ft maximum altitude above sea level |

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